

Recombinant Envelope glycoprotein protein (Cis-Tag) (STJP000147)

STJP000147

GENERAL INFORMATION

Product Type	Proteins
Short Description	Recombinant Envelope glycoprotein-Cis-Tag protein was developed from and has a target region of Cis-Tag. For use in research applications.
Host/Source	

PRODUCT PROPERTIES

Concentration	1 mg/ml
Formulation	PBS pH 7.4
Purification	His-tag, added to the C-terminal part of GP2 is used for purification of GP1-GP2 complex. Purified by metal-affinity chromatography and gel-filtration from cell culture growth media
Dilution Range	
Storage Instruction	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	
Sequence	