


## Human MT(Melatonin) ELISA Kit



### Description

The kit is a competitive enzyme immunoassay for in vitro quantitative measurement of MT in human serum, plasma and other biological fluids

<b>Model</b>	STJ150149
<b>Applications</b>	ELISA
<b>Dilution range</b>	15.625--1000pg/mL
<b>Tissue Specificity</b>	Serum, Plasma, Biological Fluids
<b>Note</b>	For Research Use Only (RUO).
<b>Molecular Weight</b>	67.4 kDa
<b>Conjugation</b>	Unconjugated
<b>Concentration</b>	1000pg/mL
<b>Storage Instruction</b>	Store at 4°C for 6 months.
<b>Function</b>	Does not bind melatonin.